



CETA Application Guide (CAG) Abstracts

as of July 19, 2022

CAG 001 - CETA Applications Guide for the use of Compounding Aseptic Isolators in Compounding Sterile Preparations in Healthcare facilities

Revised December 2008

The introduction of USP Chapter 7971 in January 2004 (revised in June 2008) along with increased pressure from some State Boards of Pharmacy has led to more robust implementation of environmental controls by facilities involved in the compounding of sterile preparations. In some cases, isolators intended for compounding sterile preparations are used either as an alternative or as an adjunct to traditional methods. Isolator use for sterile compounding is relatively new to the United States. Few isolator standards exist and none of those that are in place have been developed with pharmacy compounding in mind. This guide is intended to provide an overview of isolators and assist individuals considering their use in purchasing, installation and commissioning issues.

CAG 002 - CETA Compounding Isolator Testing Guide

Revised December 2008

The purpose of this document is to establish an industry-based minimum set of testing criteria appropriate for all Compounding Isolators used pursuant to USP Chapter 797 1. Compounding Isolators consist of Compounding Aseptic Isolators used for compounding sterile preparations and Compounding Aseptic Containment Isolators used for compounding sterile hazardous drug preparations in pharmacy applications. While this document gives general guidance or referenced guidance through relevant industry documents, it is not the intention to set the specific acceptance criteria. It is the manufacturer's responsibility to determine exact testing procedures consistent with these guidelines and assign appropriate values pertaining to acceptance criteria that is consistent with user requirements. This guide has also been established to create a uniform approach for field certifiers to allow consistent and repeatable testing at all facilities.

CAG 003 - CETA Certification Guide for Sterile Compounding Facilities

Revised May 2015

The purpose of this guide is to establish an industry-based minimum set of criteria appropriate for performance evaluation and certification of facility and environmental controls used for compounding sterile preparations. It is intended to assist compounders, facilities managers and certification



professionals in determining appropriate tests and procedures to be employed on the various engineering controls. This guideline has been established to create a uniform approach for field certifiers to allow consistent and repeatable testing at all facilities. The approach of this guide is to reference the applicable accepted industry guidelines, standards or recommended practices along with the specific tests within that document. When industry guidance documents are not available for a specific procedure, guidance will be provided here for developing an appropriate procedure.

CAG 004 - Biological Decontamination and Disinfection of Accessible Surfaces in Biosafety Cabinets

Revised September 2020

This CETA Application Guide serves as a guideline for the use of liquid chemical disinfectants to biologically disinfect (decontaminate) the readily accessible areas of Biological Safety Cabinets (BSCs). The guide serves as a general overview of these topics:

- The terms and definitions relating to surface disinfection.
- The precautions and safety measures that must be followed.
- When to surface decontaminate a BSC.
- The importance of Standard Operating Procedure (SOP) development.
- Remediation of biohazardous spills within the BSC.
- Different disinfection requirements for certifiers vs. users.
- A comparison of different microbe's susceptibility to disinfectants.
- Characteristics of an ideal disinfectant.
- Factors concerning the use of disinfectants in BSCs.
- Tradeoffs (pros and cons) of different surface disinfectants.
- Supplemental information on disinfectants.

CAG 005 - Servicing Hazardous Drug Compounding Primary Engineering Controls

Revised 2007

This document applies to hazardous drug (HD) compounding practice locations and is intended to provide guidance for service personnel such as certification technicians, facility service workers, equipment maintenance personnel, and drug compounding professionals regarding hazardous drug exposure when primary engineering controls are opened for service or repair work. This guide is based on best practice and stresses universal precautions by assuming materials are contaminated unless proven otherwise.



CAG 006 - CETA High Efficiency Filter Application Guide

Revised February 2010

The following are recommended guidelines for the selection, procurement, handling, and replacement of high efficiency filters for biosafety cabinets, isolators and laminar flow cabinets. These guidelines are intended to supplement the filter specification portion of the CETA Specification Guide for HEPA and ULPA filters procured for field replacement applications. The information enclosed is also intended to provide additional information regarding general filter guidelines for field applications.

CAG 007 - Exhaust System Requirements of Class II Biosafety Cabinets

Revised September 2020

Class II Biosafety Cabinets (BSCs) are Primary Engineering Controls (PECs), that protect personnel, product, and the environment from biohazardous aerosols and volatile toxic chemicals appropriate for use in the BSC, after performing a chemical risk assessment. If these specialized pieces of equipment are connected to building exhaust systems, there are special requirements that set them apart from chemical fume hoods and other exhausted devices. This guide is intended to help the reader understand what those special system requirements are and how a BSC should be properly connected to a building exhaust system.

CAG 008 - Certification Matrix for Sterile & Nonsterile USP Compounding Facilities

Revised March 2022

The need to have Sterile Compounding Facilities certified on a semi-annual basis has been dictated by current USP Guidelines. This Application Guide will serve as the baseline to have a Sterile Compounding Facility and its Primary Engineering Controls successfully certified in accordance with current ISO, USP, NSF and IEST standards. The requirements are laid out, test methods specified, acceptance criteria listed, and equipment needed to test described secondary engineering controls. Consideration should be taken to choose which testing to use based on Secondary and Primary Engineering Control Type.

CAG 009 - Viable Environmental Monitoring for Sterile Compounding Facilities

Revised September 2020

Viable environmental monitoring is an essential part of evaluating the state of microbial control in the sterile compounding environment. If environmental monitoring is not performed properly or frequently enough, it becomes difficult to understand the true microbial state of the tested environment. Without this understanding, patients could be at risk of receiving contaminated sterile compounded preparations. This guide provides detail on how to meet the minimum compounding standards, as well as best practice recommendations.



CAG 010 - CETA Application Guide for Informational Notes to Meet the NSF/ANSI 49:2010a Standard Requirements

Revised July 2011

The NSF/ANSI 49 Standard for Biosafety Cabinets: Design, Construction, Performance, and Field Certification had a significant number of changes and/or updates for the new 2010 release. The purpose of this document is to review the changes and provide some background and direction for compliance. In addition to the stated changes, a policy change has been implemented by NSF in the certification accreditation program to provide guidance to NSF accredited certifiers on how to apply these new standard changes. This policy change will be discussed in appendix A.

CAG 011 - Gloved Fingertip Testing for Sterile Compounding Personnel

Revised February 2021

Sterile compounding personnel pose the greatest risk to the microbial integrity of a compounded sterile preparation (CSP) through inherent bioburden and improper garbing practices. Gloved fingertip testing evaluates a compounding personnel's competency in performing hand hygiene and garbing based on the facility's standard operating procedures.

CAG 013 – Media Fill Testing for Sterile Compounding Personnel

Revised August 2021

Sterile compounding personnel pose the greatest risk to the microbial integrity of a compounded sterile preparation (CSP) through inherent bioburden and improper or poor technique when preparing and handling CSPs. Media-fill testing evaluates compounding personnel's aseptic technique and ability to safely prepare CSPs for patients.

CAG 014 – Airflow Visualization Study

Revised April 2022

Airflow visualization studies completed in sterile compounding facilities should be performed in a consistent manner and in such a way that the tests are easily repeatable and well understood by the compounding facility and certification company. CAG-014 provides guidance, methodology and industry best practices applicable to airflow visualization studies in primary engineering controls and cleanrooms.